

ScudoTex[®]



ScudoTex[®]

MEDROMA

www.medroma.co.za

ORTHOTICS

PROSTHETICS

COMPRESSION GARMENTS

GENERAL CATALOGUE



Made in Italy

IMPORTANTE NOVITÀ SCUDOTEX: LA CINTURA MULTI-GRIPPANTE REGOLABILE IN QUALSIASI PUNTO

Scudotex has recently introduced the new post-surgical, post-partum and thoracic Multi-Grippante belt, in innovative breathable hygienic fabric, with a special external surface that allows the velcro fastening to perfectly adhere at any point of the belt, offering many advantages.

Scudotex Multi-Grippante is actually particularly recommended in all cases in which an abdominal reduction occurs, such as during the post-partum or the post-surgical course of specific surgeries, offering an equilibrated back-lumbar support and an efficient pelvic-abdominal compressive action.

Another advantage of the Scudotex Multi-Grippante belt is the one to assure a perfect fit, guaranteeing at the same time the best comfort. Scudotex Multi-Grippante belt is available in the Classic or Double Closure model and it can satisfy any need in only 3 sizes, making it easy to choose the right size even when it is not possible to measure the waist circumference of the patient in bed or in hospital. Fabric and knitting 100% made in Italy.

165 BIANCA cm 28 (I-III)
163 BIANCA cm 25 (I-III)

CINTURA MULTI-GRIPPANTE REGOLABILE DOPPIA CHIUSURA

Breathable elastic belts made in elasticized fabric, with a special external surface knitting that allows the velcro double closure to perfectly adhere at any point of the belt, offering the best fit and comfort.

CON DOPPIA CHIUSURA

*in only
3 sizes*


155 BIANCA cm 28 (I-III)
153 BIANCA cm 25 (I-III)
145 BIANCA cm 20 (I-III)
143 BIANCA cm 16 (I-III)

CINTURA MULTI-GRIPPANTE REGOLABILE CLASSICA

Breathable elastic belts made in elasticized fabric, with a special external surface knitting that allows the velcro closure to perfectly adhere at any point of the belt, offering the best fit and comfort. The heights cm 20 and cm 16 are also recommended as thoracic belt.

Misura	I	II	III
Size	S/M	L/XL	2XL/3XL
Circonferenza vita / torace (cm)	60-90	80-120	110-150
Waist / chest measurement (cm)			

MEDIUM COMPRESSION mm Hg 18–21



451 GAMBALETTO UOMO SILVER SKY COTTON

Graduated compression elastic knee-high, in hygienic and fresh cotton, opaque, with bacteriostatic microfiber with silver ions that guarantees a high transpiration and comfort, maintaining the natural bacterial balance of the skin. Elegant ribbed knitting, particularly appreciated, also in the summer season, by people who carry out sedentary activities or stand many hours in very small spaces.

Available in colours:
antracite, blu and nero (II-VI).

CINTURE SPECIALI

CINTURA PER COLOSTOMIA

pannello staccabile, foro Ø cm 7 e 10

Special elastic belt in cotton which offers a balanced lumbosacral and abdominal support, recommended to colostomates.

The abdominal padded panel, fixed on the elastic belt, has a waterproof surface against the skin for a hygienic cleaning and a semirigid ring for the fastening of the pouch to the stoma.

The abdominal panel can be bought also separately.

To be worn according to physician's advice.



**C761 CINTURA COLOSTOMIA cm 24,
CON PANNELLO Ø cm 7 (I-VII)**

C760 CINTURA COLOSTOMIA cm 24,

Misura - Size	I (XS)	II (S)	III (M)	IV (L)	V (XL)	VI (XXL)	VII (XXXL)
Circonf. vita cm Waist circum. cm	60 - 75	70 - 85	80 - 95	90 - 105	100 - 115	110 - 125	120 - 135



CORSETTI STECCATI TUTORI ORTOPEDICI LOMBOSACRALI

Orthopedic stayed braces are very important in the prevention and care of some spinal column pathologies, with particular reference to the lumbosacral area mainly stressed by external strains following physical and lifting activities. Thanks to the regular and supervised use of the orthopedic corset, the rachis motility can be reduced and a better balance of the vertebral weight can be obtained, with positive effects on the posture. The use of the orthopedic corsets is recommended to people suffering from low back pain, sciatica, discopathy, or generic backache and to people who subject their rachis to frequent efforts. From this page on you will find a complete orthopedic corsets line so that each patient can choose the most appropriate model according to the physician's indications. All the lumbosacral corsets shown in these pages are recommended for work, sport and to avoid labour accidents, for the care of the worker's health



CORSETTO ORTOPEDICO SUPER CROSS

The Super Cross Corset is an orthopedic brace in exclusive hygienic perspiring elastic cotton fabric, reinforced by a back crossed support and a double abdominal band that shapes the body and perfectly fits both man and woman. Designed with 4 semirigid removable paravertebral stays and abdominal antirolling stays with adjustable adhesive closure. Available in the following colors and heights.



610 NATURALE cm 27 (I-VIII)
618 NATURALE cm 24 (II-VI)
667 NATURALE cm 32 (II-VI)



Misura	I	II	III	IV	V	VI	VII	VIII
Size	XS	S	M	L	XL	XXL	XXXL	XXXXL
Circonferenza vita (cm) Waist circumference (cm)	60-72	70-92	90-104	100-116	110-128	120-140	130-152	140-165

TUTORI ORTOPEDICI STECCATI



GINOCCHIERA ARTICOLATA STABILIZZANTE 844 cm 30 - 846 cm 40

The Rekordsan hinged knee brace is an orthopedic brace with open knee cap area, made in "high density" neoprene fabric that accumulates the body temperature and warms efficiently the articulation. Thanks to its supporting and containing action, the blood circulation is improved, muscular contractions and swelling are reduced. It has four monocentric metal stays.

To be worn in accordance with medical advise in cases of: • dislocations • distortions • luxations • light or moderated ligaments instabilities • tendonitis • ligaments inflammation • rheumatic illnesses • arthritis • arthrosis • after surgical operation • degenerative changes.

Misura Size	I	II	III	IV
	S	M	L	XL
Circonferenza sopra ginocchio (cm) Above knee measurement (cm)	30-34	34-38	38-42	42-46

CALZE ELASTICHE MAGLIA A RETE A COMPRESSIONE GRADUATA TERAPEUTICA

140 DENARI MAGLIA A RETE COMPRESSIONE MEDIO-FORTE mm Hg 19-22

Characteristics: Scudotex elastic hosiery manufactured with very high quality yarns, net knitting, perspiring, reinforced heel, flat seams. Comfortable elastic hosiery, well tolerated, develops a healthy massaging and toning up action on the legs and assures good results in the treatment of the Chronic Venous Insufficiency.

Compression and indications: graduated medium-strong compression hosiery, mm Hg 19-22, curative, recommended for the treatment of the venous insufficiency, also chronic, with possible development of varicose veins and formation of oedema. Elegant, lightly opaque, it hides the slight imperfection of the skin. To be worn according to physician's advice.

- 490 COLLANT 140 DENARI EXTRA
- 493 COLLANT 140 DENARI CLASSIC
- 496 GESTANTE 140 DENARI EXTRA
- 489 AUTOREGGENTE 140 DENARI EXTRA
- 497 CALZA 140 DENARI EXTRA
- 498 GAMBALETTA 140 DENARI EXTRA

anatomic panty, knitted waistband, with a big gusset from the 2nd size (I-VI).
 anatomic panty, applied waistband, with a big gusset from the 2nd size (II-VI).
 pantyhose for pregnancy with a large polyextensible front, for an equilibrated abdominal support (II- VI).
 thigh high with an attractive lace and siliconated anti-slipping edge (I-VI).
 stocking model with a special soft edge, for an easy fastening of the suspender belt (I-VI).
 knee high with a special soft elastic (I-VI).